





COLLECTION KIT	CONTAINER	TEST NAMES
<p><b>Aptima® Multitest Swab</b></p> <p>Store specimen at <u>room temperature</u> after collection.</p>		<p>Aptima® Combo 2 Assay Chlamydia trachomatis &amp; Neisseria gonorrhoeae</p> <p>Aptima® Mycoplasma genitalium Assay</p> <p>Aptima® Trichomonas vaginalis Assay</p> <p>Aptima® Bacterial vaginosis Assay</p> <p>Aptima® Candida vaginosis &amp; Trichomonas vaginalis Assay</p> <p>Aptima® Herpes Simplex Virus, 1 &amp; 2</p>
<p><b>Aptima® Urine Specimen</b></p> <p>Store specimen at <u>room temperature</u> after collection.</p> <p>Patient should not urinate for at least 1 hour prior to specimen collection.</p> <p>Send only the Aptima® Urine Specimen Tube. <u>No urine cups will be accepted.</u></p> <ul style="list-style-type: none"> <li>• Direct patient to provide first-catch urine (~20-30 mL of initial urine stream) into urine collection cup. Collection of larger volumes of urine may result in specimen dilution that may reduce test sensitivity. Female patients should not cleanse labial area prior to providing specimen.</li> <li>• Remove cap and transfer 2 mL of urine into urine specimen transport tube using the pipette provided.</li> <li>• The correct volume of urine has been added when the fluid level is between black fill lines on urine specimen transport tube label.</li> <li>• Discard pipette. Tightly screw cap onto tube. This is now known as the "processed urine specimen."</li> </ul>		<p><b><u>Female Urine</u></b></p> <p>Aptima® Combo 2 Assay Chlamydia trachomatis &amp; Neisseria gonorrhoeae</p> <p>Aptima® Mycoplasma genitalium Assay</p> <p>Aptima® Trichomonas vaginalis Assay</p> <p><b><u>Male Urine</u></b></p> <p>Aptima® Combo 2 Assay Chlamydia trachomatis &amp; Neisseria gonorrhoeae</p> <p>Aptima® Mycoplasma genitalium Assay</p> <p>Aptima® Trichomonas vaginalis Assay</p>
<p><b>Aptima® Unisex Swab</b></p> <p>Store specimen at <u>room temperature</u> after collection.</p> <p>Patient should not urinate for at least 1 hour prior to specimen collection.</p> <p><b><u>Male Urethral Swab Specimens</u></b></p> <ul style="list-style-type: none"> <li>• Discard white cleaning swab. It is not needed for male specimen collection.</li> <li>• Insert blue specimen collection swab 2-4 cm into urethra. Gently rotate swab clockwise for 2-3 seconds. Withdraw swab. If the soft tip is touched, laid down, or dropped, discard and replace with a new Aptima Unisex Swab Specimen Collection Kit.</li> <li>• Immediately place swab into transport tube so black score line is at top of tube. Align the score line with top edge of tube and carefully break shaft. Discard top portion of shaft. Tightly screw cap onto tube.</li> </ul> <p><b><u>Female Endocervical Swab Specimens</u></b></p> <ul style="list-style-type: none"> <li>• Use white cleaning swab to remove excess mucus from cervical os and surrounding mucosa. Discard.</li> <li>• Insert blue collection swab into endocervical canal. Gently rotate swab clockwise for 10-30 seconds. Withdraw swab; avoid any contact with vaginal mucosa. If the soft tip is touched, laid down, or dropped, discard and get a new Aptima Unisex Swab Specimen Collection Kit.</li> <li>• While holding swab in hand, unscrew tube cap. Do not spill tube contents. If spilled, discard and replace with a new Aptima Unisex Swab Specimen Kit.</li> <li>• Immediately place swab into transport tube so black score line is at top of tube. Align the score line with top edge of tube and carefully break shaft. Discard top portion of shaft. Tightly screw cap onto tube.</li> </ul>		<p>Aptima® Combo 2 Assay Chlamydia trachomatis &amp; Neisseria gonorrhoeae</p> <p>Aptima® Mycoplasma genitalium Assay</p> <p>Aptima® Trichomonas vaginalis Assay</p>
<p><b>ThinPrep® Pap Test Vial</b></p> <p>Store specimen at <u>room temperature</u> after collection.</p>		<p>Aptima® Combo 2 Assay Chlamydia trachomatis &amp; Neisseria gonorrhoeae</p> <p>Aptima® Trichomonas vaginalis Assay</p> <p>Aptima® HPV Assay</p> <p>Aptima® HPV 16 18/45 Genotype Assay</p>