

Urgent Medical Device Recall • BD Microtainer® Tubes

Notification Date: 17 MAY 2019

Effective Date: 17 MAY 2019

Urgent Medical Device Recall – BD Microtainer® Tubes

The South Bend Medical Foundation has been notified that BD is conducting a voluntary medical device recall for various lots of BD Microtainer® Tubes based on confirmation that there may be damaged tube reservoirs present. A damaged reservoir may lead to a decreased fill volume causing samples to be insufficient for testing and an improper blood-to-additive ratio, potentially producing erroneous results. Additionally, a damaged tube reservoir may result in the potential for blood leakage and subsequent exposure to healthcare workers. Distribution of the affected lots began on **April 11, 2018**.

The root cause investigation has determined that the damage to the tube reservoirs was caused by an equipment malfunction during the molding process. The equipment malfunction has been corrected to prevent recurrence of the identified root cause.

Immediately review your inventory for the specific catalog and lot numbers listed below. Please return all affected tubes back to **STOCKROOM, CLAB** no later than **Monday, May 27, 2019** and mark your return as **BD RECALL**.

Please utilize our online supply order form on the sbmf.org Client Resources page when ordering replacement tubes. You may also e-mail Stockroom directly at Stockroom@sbmf.org. If you need help placing your new supply order, please contact our Client Services Department at 800-950-7263.

Please use the BD contact information below for complaints, adverse event reports or questions regarding this recall

BD Contact Department	US Contact Information
Customer Quality	888-237-2762 OPT 3, OPT 2 Monday – Friday 8:00am and 5:00pm (CT)

[SBMF online Test Directory](#)

Questions: Please contact **CLIENT SERVICES 800-950-7263**

Technical Notices and Lab Alerts are distributed electronically.

Email addresses may be added/unsubscribed at our website: [Client Notices](#)

SOUTH BEND MEDICAL FOUNDATION

530 North Lafayette Boulevard • South Bend, Indiana 46601
(574) 234-4176 • (800) 544-0925 • www.sbmf.org
Joyce L. Simpson, M.D. • Medical Director

Product	Part / Catalog Number	Lot / Serial Number	UDI (GTIN, DI + PI)	Exp. Date
BD Microtainer® Z (No Additive Tubes)	365963	8117750	(01)30382903659631 (17)191031(10)8117750(30)50 (17)191031(10)8117750(30)200 (01)50382903659635	2019-10-31
BD Microtainer® Z (No Additive Tubes)	365963	8155792	(01)30382903659631 (17)191130(10)8155792(30)50 (17)191130(10)8155792(30)200 (01)50382903659635	2019-11-30
BD Microtainer® Z (No Additive Tubes)	365963	8223778	(01)30382903659631 (17)200131(10)8223778(30)50 (17)200131(10)8223778(30)200 (01)50382903659635	2020-01-31
BD Microtainer® Z (No Additive Tubes)	365963	8223779	(01)30382903659631 (17)200131(10)8223779(30)50 (17)200131(10)8223779(30)200 (01)50382903659635	2020-01-31
BD Microtainer® Z (No Additive Tubes)	365963	8268514	(01)30382903659631 (17)200229(10)8268514(30)50 (17)200229(10)8268514(30)200 (01)50382903659635	2020-02-29
BD Microtainer® Tubes with LH (Lithium Heparin)	365965	811767N	(01)30382903659655 (17)190930(10)811767N(30)50 (17)190930(10)811767N(30)200 (01)50382903659659	2019-09-30
BD Microtainer® Tubes with LH (Lithium Heparin)	365965	815066N	(01)30382903659655 (17)191130(10)815066N(30)50 (17)191130(10)815066N(30)200 (01)50382903659659	2019-11-30
BD Microtainer® Tubes with LH (Lithium Heparin)	365965	821471N	(01)30382903659655 (17)200131(10)821471N(30)50 (17)200131(10)821471N(30)200 (01)50382903659659	2020-01-31
BD Microtainer® Tubes with LH (Lithium Heparin)	365965	821990N	(01)30382903659655 (17)200131(10)821990N(30)50 (17)200131(10)821990N(30)200 (01)50382903659659	2020-01-31
BD Microtainer® SST™	365967	8094573	(01)30382903659679 (17)190630(10)8094573(30)50 (17)190630(10)8094573(30)200 (01)50382903659673	2019-06-30
BD Microtainer® SST™	365967	8094575	(01)30382903659679 (17)190630(10)8094575(30)50 (17)190630(10)8094575(30)200 (01)50382903659673	2019-06-30
BD Microtainer® SST™	365967	8094578	(01)30382903659679 (17)190630(10)8094578(30)50 (17)190630(10)8094578(30)200 (01)50382903659673	2019-06-30
BD Microtainer® SST™	365967	8117697	(01)30382903659679 (17)190731(10)8117697(30)50 (17)190731(10)8117697(30)200 (01)50382903659673	2019-07-31
BD Microtainer® SST™	365967	8165883	(01)30382903659679 (17)190831(10)8165883(30)50 (17)190831(10)8165883(30)200 (01)50382903659673	2019-08-31
BD Microtainer® SST™	365967	8165884	(01)30382903659679 (17)190831(10)8165884(30)50 (17)190831(10)8165884(30)200 (01)50382903659673	2019-08-31
BD Microtainer® SST™	365967	8170911	(01)30382903659679 (17)190831(10)8170911(30)50 (17)190831(10)8170911(30)200 (01)50382903659673	2019-08-31
BD Microtainer® SST™	365967	8171861	(01)30382903659679 (17)190831(10)8171861(30)50 (17)190831(10)8171861(30)200 (01)50382903659673	2019-08-31
BD Microtainer® SST™	365967	8220937	(01)30382903659679 (17)191031(10)8220937(30)50 (17)191031(10)8220937(30)200 (01)50382903659673	2019-10-31
BD Microtainer® SST™	365967	8220938	(01)30382903659679 (17)191031(10)8220938(30)50 (17)191031(10)8220938(30)200 (01)50382903659673	2019-10-31
BD Microtainer® SST™	365967	8282848	(01)30382903659679 (17)191231(10)8282848(30)50 (17)191231(10)8282848(30)200 (01)50382903659673	2019-12-31
BD Microtainer® SST™	365967	8284571	(01)30382903659679 (17)191231(10)8284571(30)50 (17)191231(10)8284571(30)200 (01)50382903659673	2019-12-31
BD Microtainer® SST™	365967	8284574	(01)30382903659679 (17)191231(10)8284574(30)50 (17)191231(10)8284574(30)200 (01)50382903659673	2019-12-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8101695	(01)30382903659747 (17)190930(10)8081695(30)50 (17)190930(10)8101695(30)200 (01)50382903659741	2019-09-30
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8101696	(01)30382903659747 (17)190930(10)8081696(30)50 (17)190930(10)8101696(30)200 (01)50382903659741	2019-09-30
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8101699	(01)30382903659747 (17)190930(10)8101699(30)50 (17)190930(10)8101699(30)200 (01)50382903659741	2019-09-30
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8101701	(01)30382903659747 (17)190930(10)8101701(30)50 (17)190930(10)8101701(30)200 (01)50382903659741	2019-09-30
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8106835	(01)30382903659747 (17)190930(10)8106835(30)50 (17)190930(10)8106835(30)200 (01)50382903659741	2019-09-30
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8106837	(01)30382903659747 (17)190930(10)8106837(30)50 (17)190930(10)8106837(30)200 (01)50382903659741	2019-09-30
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8106842	(01)30382903659747 (17)190930(10)8106842(30)50 (17)190930(10)8106842(30)200 (01)50382903659741	2019-09-30

SOUTH BEND MEDICAL FOUNDATION

530 North Lafayette Boulevard • South Bend, Indiana 46601
 (574) 234-4176 • (800) 544-0925 • www.sbmf.org
 Joyce L. Simpson, M.D. • Medical Director

Product	Part / Catalog Number	Lot / Serial Number	UDI (GTIN, DI + PI)	Exp. Date
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8106847	(01)30382903659747 (17)190930(10)8106847(30)50 (17)190930(10)8106847(30)200 (01)50382903659741	2019-09-30
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8106850	(01)30382903659747 (17)190930(10)8106850(30)50 (17)190930(10)8106850(30)200 (01)50382903659741	2019-09-30
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8120850	(01)30382903659747 (17)191031(10)8120850(30)50 (17)191031(10)8120850(30)200 (01)50382903659741	2019-10-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8120851	(01)30382903659747 (17)191031(10)8120851(30)50 (17)191031(10)8120851(30)200 (01)50382903659741	2019-10-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8121979	(01)30382903659747 (17)191031(10)8121979(30)50 (17)191031(10)8121979(30)200 (01)50382903659741	2019-10-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8121981	(01)30382903659747 (17)191031(10)8121981(30)50 (17)191031(10)8121981(30)200 (01)50382903659741	2019-10-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8124783	(01)30382903659747 (17)191031(10)8124783(30)50 (17)191031(10)8124783(30)200 (01)50382903659741	2019-10-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8124880	(01)30382903659747 (17)191031(10)8124880(30)50 (17)191031(10)8124880(30)200 (01)50382903659741	2019-10-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8127770	(01)30382903659747 (17)191031(10)8127770(30)50 (17)191031(10)8127770(30)200 (01)50382903659741	2019-10-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8128575	(01)30382903659747 (17)191031(10)8128575(30)50 (17)191031(10)8128575(30)200 (01)50382903659741	2019-10-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8129642	(01)30382903659747 (17)191031(10)8129642(30)50 (17)191031(10)8129642(30)200 (01)50382903659741	2019-10-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8129645	(01)30382903659747 (17)191031(10)8129645(30)50 (17)191031(10)8129645(30)200 (01)50382903659741	2019-10-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8197705	(01)30382903659747 (17)191231(10)8197705(30)50 (17)191231(10)8197705(30)200 (01)50382903659741	2019-12-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8197970	(01)30382903659747 (17)191231(10)8197970(30)50 (17)191231(10)8197970(30)200 (01)50382903659741	2019-12-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8198921	(01)30382903659747 (17)191231(10)8198921(30)50 (17)191231(10)8198921(30)200 (01)50382903659741	2019-12-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8198927	(01)30382903659747 (17)191231(10)8198927(30)50 (17)191231(10)8198927(30)200 (01)50382903659741	2019-12-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8243545	(01)30382903659747 (17)200229(10)8243545(30)50 (17)200229(10)8243545(30)200 (01)50382903659741	2020-02-29
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8247628	(01)30382903659747 (17)200229(10)8247628(30)50 (17)200229(10)8247628(30)200 (01)50382903659741	2020-02-29
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8248705	(01)30382903659747 (17)200229(10)8248705(30)50 (17)200229(10)8248705(30)200 (01)50382903659741	2020-02-29
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8249986	(01)30382903659747 (17)200229(10)8249986(30)50 (17)200229(10)8249986(30)200 (01)50382903659741	2020-02-29
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8288904	(01)30382903659747 (17)200331(10)8288904(30)50 (17)200331(10)8288904(30)200 (01)50382903659741	2020-03-31
BD Microtainer® Tubes with K2E (K2EDTA)	365974	8290928	(01)30382903659747 (17)200331(10)8290928(30)50 (17)200331(10)8290928(30)200 (01)50382903659741	2020-03-31
BD Microtainer® SST™ - Amber	365978	8193502	(01)30382903659785 (17)190930(10)8193502(30)50 (17)190930(10)8193502(30)200 (01)50382903659789	2019-09-30

SOUTH BEND MEDICAL FOUNDATION

530 North Lafayette Boulevard • South Bend, Indiana 46601
(574) 234-4176 • (800) 544-0925 • www.sbmf.org
Joyce L. Simpson, M.D. • Medical Director

Product	Part / Catalog Number	Lot / Serial Number	UDI (GTIN, DI + PI)	Exp. Date
BD Microtainer® SST™ - Amber	365978	8193507	(01)30382903659785 (17)190930(10)8193507(30)50 (17)190930(10)8193507(30)200 (01)50382903659789	2019-09-30
BD Microtainer® SST™ - Amber	365978	8193511	(01)30382903659785 (17)190930(10)8193511(30)50 (17)190930(10)8193511(30)200 (01)50382903659789	2019-09-30
BD Microtainer® SST™ - Amber	365978	8204927	(01)30382903659785 (17)190930(10)8204927(30)50 (17)190930(10)8204927(30)200 (01)50382903659789	2019-09-30
BD Microtainer® SST™ - Amber	365978	8205758	(01)30382903659785 (17)190930(10)8205758(30)50 (17)190930(10)8205758(30)200 (01)50382903659789	2019-09-30
BD Microtainer® SST™ - Amber	365978	8207774	(01)30382903659785 (17)190930(10)8207774(30)50 (17)190930(10)8207774(30)200 (01)50382903659789	2019-09-30
BD Microtainer® SST™ - Amber	365978	8208542	(01)30382903659785 (17)190930(10)8208542(30)50 (17)190930(10)8208542(30)200 (01)50382903659789	2019-09-30
BD Microtainer® SST™ - Amber	365978	8208544	(01)30382903659785 (17)190930(10)8208544(30)50 (17)190930(10)8208544(30)200 (01)50382903659789	2019-09-30
BD Microtainer® SST™ - Amber	365978	8232513	(01)30382903659785 (17)191031(10)8232513(30)50 (17)191031(10)8232513(30)200 (01)50382903659789	2019-10-31
BD Microtainer® SST™ - Amber	365978	8233967	(01)30382903659785 (17)191031(10)8233967(30)50 (17)191031(10)8233967(30)200 (01)50382903659789	2019-10-31
BD Microtainer® SST™ - Amber	365978	8236755	(01)30382903659785 (17)191031(10)8236755(30)50 (17)191031(10)8236755(30)200 (01)50382903659789	2019-10-31
BD Microtainer® SST™ - Amber	365978	8236756	(01)30382903659785 (17)191031(10)8236756(30)50 (17)191031(10)8233967(30)200 (01)50382903659789	2019-10-31
BD Microtainer® SST™ - Amber	365978	8239857	(01)30382903659785 (17)191031(10)8239857(30)50 (17)191031(10)8239857(30)200 (01)50382903659789	2019-10-31
BD Microtainer® SST™ - Amber	365978	8260948	(01)30382903659785 (17)191130(10)8260948(30)50 (17)191130(10)8260948(30)200 (01)50382903659789	2019-11-30
BD Microtainer® SST™ - Amber	365978	8262866	(01)30382903659785 (17)191130(10)8262866(30)50 (17)191130(10)8262866(30)200 (01)50382903659789	2019-11-30
BD Microtainer® SST™ - Amber	365978	8282812	(01)30382903659792 (17)191231(10)8282812(30)50 (17)191231(10)8282812(30)200 (01)50382903659796	2019-12-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	808592N	(01)30382903659853 (17)190630(10)808592N(30)50 (17)190630(10)808592N(30)200 (01)50382903659857	2019-06-30
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	809669N	(01)30382903659853 (17)190630(10)809669N(30)50 (17)190630(10)809669N(30)200 (01)50382903659857	2019-06-30
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	811070N	(01)30382903659853 (17)190630(10)811070N(30)50 (17)190630(10)811070N(30)200 (01)50382903659857	2019-06-30
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	811071N	(01)30382903659853 (17)190630(10)811071N(30)50 (17)190630(10)811071N(30)200 (01)50382903659857	2019-06-30
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	811664N	(01)30382903659853 (17)190630(10)811664N(30)50 (17)190630(10)811664N(30)200 (01)50382903659857	2019-06-30
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	813767N	(01)30382903659853 (17)190731(10)813767N(30)50 (17)190731(10)813767N(30)200 (01)50382903659857	2019-07-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	815060N	(01)30382903659853 (17)190731(10)815060N(30)50 (17)190731(10)815060N(30)200 (01)50382903659857	2019-07-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	815061N	(01)30382903659853 (17)190831(10)815061N(30)50 (17)190831(10)815061N(30)200 (01)50382903659857	2019-08-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	815062N	(01)30382903659853 (17)190831(10)815062N(30)50 (17)190831(10)815062N(30)200 (01)50382903659857	2019-08-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	816499N	(01)30382903659853 (17)190831(10)816499N(30)50 (17)190831(10)816499N(30)200 (01)50382903659857	2019-08-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	821563N	(01)30382903659853 (17)191031(10)821563N(30)50 (17)191031(10)821563N(30)200 (01)50382903659857	2019-10-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	821564N	(01)30382903659853 (17)191031(10)821564N(30)50 (17)191031(10)821564N(30)200 (01)50382903659857	2019-10-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	821565N	(01)30382903659853 (17)191031(10)821565N(30)50 (17)191031(10)821565N(30)200 (01)50382903659857	2019-10-31

SOUTH BEND MEDICAL FOUNDATION

530 North Lafayette Boulevard • South Bend, Indiana 46601
 (574) 234-4176 • (800) 544-0925 • www.sbmf.org
 Joyce L. Simpson, M.D. • Medical Director

Product	Part / Catalog Number	Lot / Serial Number	UDI (GTIN, DI + PI)	Exp. Date
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	821888N	(01)30382903659853 (17)191031(10)821888N(30)50 (17)191031(10)821888N(30)200 (01)50382903659857	2019-10-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	827063N	(01)30382903659853 (17)191231(10)827063N(30)50 (17)191231(10)827063N(30)200 (01)50382903659857	2019-12-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	827064N	(01)30382903659853 (17)191231(10)827064N(30)50 (17)191231(10)827064N(30)200 (01)50382903659857	2019-12-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	827589N	(01)30382903659853 (17)191231(10)827589N(30)50 (17)191231(10)827589N(30)200 (01)50382903659857	2019-12-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	827693N	(01)30382903659853 (17)191231(10)827693N(30)50 (17)191231(10)827693N(30)200 (01)50382903659857	2019-12-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin)	365985	827699N	(01)30382903659853 (17)191231(10)827699N(30)50 (17)191231(10)827699N(30)200 (01)50382903659857	2019-12-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	819352N	(01)30382903659877 (17)190930(10)819352N(30)50 (17)190930(10)819352N(30)200 (01)50382903659871	2019-09-30
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	822851N	(01)30382903659877 (17)191031(10)822851N(30)50 (17)191031(10)822851N(30)200 (01)50382903659871	2019-10-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	822852N	(01)30382903659877 (17)191031(10)822852N(30)50 (17)191031(10)822852N(30)200 (01)50382903659871	2019-10-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	822853N	(01)30382903659877 (17)191031(10)822853N(30)50 (17)191031(10)822853N(30)200 (01)50382903659871	2019-10-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	822854N	(01)30382903659877 (17)191031(10)822854N(30)50 (17)191031(10)822854N(30)200 (01)50382903659871	2019-10-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	822855N	(01)30382903659877 (17)191031(10)822855N(30)50 (17)191031(10)822855N(30)200 (01)50382903659871	2019-10-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	822862N	(01)30382903659877 (17)191031(10)822862N(30)50 (17)191031(10)822862N(30)200 (01)50382903659871	2019-10-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	824092N	(01)30382903659877 (17)191031(10)824092N(30)50 (17)191031(10)824092N(30)200 (01)50382903659871	2019-10-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	826768N	(01)30382903659877 (17)191130(10)826768N(30)50 (17)191130(10)826768N(30)200 (01)50382903659871	2019-11-30
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	827652N	(01)30382903659877 (17)191231(10)827652N(30)50 (17)191231(10)827652N(30)200 (01)50382903659871	2019-12-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	827753N	(01)30382903659877 (17)191231(10)827753N(30)50 (17)191231(10)827753N(30)200 (01)50382903659871	2019-12-31
BD Microtainer® PST™ Tubes with LH (Lithium Heparin) - Amber	365987	827754N	(01)30382903659877 (17)191231(10)827754N(30)50 (17)191231(10)827754N(30)200 (01)50382903659871	2019-12-31

SOUTH BEND MEDICAL FOUNDATION

530 North Lafayette Boulevard • South Bend, Indiana 46601
 (574) 234-4176 • (800) 544-0925 • www.sbmf.org
 Joyce L. Simpson, M.D. • Medical Director

Product	Part / Catalog Number	Lot / Serial Number	UDI (GTIN, DI + PI)	Exp. Date
BD Microtainer® Tubes with FE (Sodium Fluoride / Disodium EDTA) Glycolytic Inhibitor	365992	8212777	(01)30382903659921 (17)191031(10)8212777(30)50 (17)191031(10)8212777(30)200 (01)50382903659925	2019-10-31
BD Microtainer® Tubes with FE (Sodium Fluoride / Disodium EDTA) Glycolytic Inhibitor	365992	8212785	(01)30382903659921 (17)191031(10)8212785(30)50 (17)191031(10)8212785(30)200 (01)50382903659925	2019-10-31
BD Microtainer® Tubes with FE (Sodium Fluoride / Disodium EDTA) Glycolytic Inhibitor	365992	8268546	(01)30382903659921 (17)191130(10)8268546(30)50 (17)191031(10)8268546(30)200 (01)50382903659925	2019-11-30
BD Microtainer® Tubes with FE (Sodium Fluoride / Disodium EDTA) Glycolytic Inhibitor	365992	8268549	(01)30382903659921 (17)191130(10)8268549(30)50 (17)191130(10)8268549(30)200 (01)50382903659925	2019-11-30

SOUTH BEND MEDICAL FOUNDATION

530 North Lafayette Boulevard • South Bend, Indiana 46601
 (574) 234-4176 • (800) 544-0925 • www.sbmf.org
 Joyce L. Simpson, M.D. • Medical Director